

SHIPPING DETAILS:

PhoThera 600

Model: PhoThera 600 with 4 or 6 Lamps



Weight:
80 lbs



Shipping Method:
Standard Ground



Container:
Cardboard Box




Unpacking Hints:
If possible, use a dolly to move item to an area large enough to accommodate device. Sit device on soft surface to avoid scratches. Remove staples to open box; caution, staples are very sharp.



Helpful Tip:
Check off items on this sheet as you empty contents.

Dimensions:

9" W x 77" L x 25" H



LIT-00063 [0]



(216)-831-0600
Press 1 for Daavlin or 2 for NatBio



info@phoThera.com