

SHIPPING DETAILS:

Phothera 400



Weight:

10 Lamp = 41 lbs
20 Lamp = 48 lbs



Shipping Method:

UPS Ground or FedEx



Container:

Cardboard Box



Unpacking Hints:

Carefully cut packaging tape to avoid damaging interior contents. The feet to the unit are contained within the packaging; please check foam before discarding.




Helpful Tip:

Your 400 device can be mounted multiple ways; either by the feet that come with the device, or by using a VESA tv wall mount suitable for the weight of the device. Check off items on this sheet as you empty contents.

LIT-00065 [0]

Dimensions:
26" W x 55" L x 9" H



Packaging Foam

Contents:

- Power Cable
- Safety Goggles

Folder Contents:

- Owner's Manual
- Warranty Information
- Time Chart
- Calibration Sheet

400 Device Under Protective Layer

Packaging Foam



(216)-831-0600

Press 1 for Daavlin or 2 for NatBio



info@phothera.com