

Phothera

100



Phothera 100

Compact, effective, and affordable. The Phothera 100 was designed to effectively treat the scalp and small areas while still maintaining a size that is compact, travels well, and is easy to store when not in use.

Phothera 100



- 1 Thin, lightweight, and compact, the Phothera 100 is an ideal option for those who need an effective device that travels and stores well.
- 2 The onboard DT Controller makes treatment entry quick and easy with a control panel that is intuitively designed, while the backlit LCD screen clearly displays the entered treatment time.
- 3 The high-intensity lamp clears as effectively as the powerful, top-of-the-line products found in clinics.
 - The Phothera 100 can be equipped with a Narrowband UVB or UVA lamp.
- 4 The durable slide-on comb attachment separates the hair for effective scalp penetration and doubles as a distance to maintain proper height from the treatment area - ensuring consistent treatments.
- 5 A nylon case conveniently stores the device and accessories in one central location, making travel and storage quick and easy.

Treatment Area	1" x 4.5"	2.54 cm x 11.4 cm
Wand Size	8.25" x 1.75"	21 cm x 5 cm
Wand Weight	4 oz	0.1 kg
Fixture Size	9" x 5"	23 cm x 13 cm
Total Height	2.75"	7 cm
Total Weight	1 lb 12 oz	0.8 kg

For further information, please contact Phothera:
216.831.0600
phothera.com | info@phothera.com