

# Phoothera 400



## Phoothera 400

The Phoothera 400 is a mid-size phototherapy device used to treat localized to widespread conditions. The light weight and portability of the device make it perfect for small rooms, apartments, and dorms. The device can be configured with variable lamp quantities and types.

Phoothera



# Phothera 400



- 1 The Phothera 400 is a robust home phototherapy device featuring a durable, yet lightweight design. The smooth acrylic front panel that protects the user from contact with the lamps and is easy to clean.
- 2 The slim lines of the Phothera 400 make it ideal for smaller rooms where space is a concern. When not in use, the device can easily be moved into a closet for storage.
- 3 High output PL-L lamps provide the power of two lamps in the space of one – delivering surprisingly short treatment times. The Phothera 400 offers choice in lamp quantity and wavelength:
  - 10 or 20 lamp
  - Narrowband UVB, UVA, or UVA1
- 4 Two injection molded platforms ensure stability and safety across a variety of common household surfaces.
- 5 The ClearLink Controller provides patient safety, security, choice, and provider peace-of-mind, packed in an intuitive touchscreen display.

Treatment Area	18" x 46.25"	45.7cm x 117.5 cm
Width (including stand)	21.5"	54.6 cm
Height	49.5"	125.7 cm
Depth	(without stand) 3.5" (with stand) 22"	(without stand) 8.9 cm (with stand) 55.9 cm
Weight	(10 Lamp) 30 lbs (20 Lamp) 37 lbs	(10 Lamp) 13.6 kgs (20 Lamp) 16.8 kgs

For further information, please contact Phothera:  
216.831.0600  
phothera.com | info@phothera.com