

Phothera 600



Phothera 600

The Phothera 600 is a full-body phototherapy device used to treat widespread conditions. The device offers a slim design and minimal footprint, is equipped with four or six lamps, and is available with either the DT or ClearLink Controller.

Phothera

Phothera 600



- 1 The Phothera 600 offers patients and prescribers choices when choosing lamp quantities and wavelengths:
 - 4 or 6 lamps
 - Narrowband UVB
- 2 The slim lines of the Phothera 600 make it ideal for a home setting; easily fitting within just 3 square feet of floor space.
- 3 Four (4) or six (6) high output, six foot lamps deliver effective phototherapeutic relief.
- 4 The Phothera 600 offers both simplicity and flexibility, with two customizable controller options designed to fit each patient's unique needs and goals.
 - ClearLink Controller - Safety, security, and patient monitoring, all in an intuitive, touch-screen controller.
 - DT Controller - Simplicity in the form of a time-only controller.

Treatment Area	19" x 70.25"	48.3 cm x 178.4 cm
Width	20"	50.8 cm
Height	74.5"	189.2 cm
Depth (including feet)	24"	61 cm
Weight	75 lbs	34 kg

For further information, please contact Phothera:
216.831.0600
phothera.com | info@phothera.com