


Phothera

100XL



Phothera 100XL

The Phothera 100XL is a folding, handheld phototherapy device that provides fast and effective targeted treatments for psoriasis, vitiligo, atopic dermatitis, and other photoresponsive conditions.

Phothera 100XL



- 1 The Phothera 100XL is a thin and lightweight handheld device featuring a collapsible design, making storing or traveling with the device incredibly convenient.
- 2 The Phothera 100XL is equipped with two Narrowband UVB lamps. These high-intensity lamps treat as effectively as the powerful, top-of-the line products found in clinics to help patients achieve clear skin.
- 3 The onboard DT Controller makes treatment entry quick and easy with a control panel that is intuitively designed, while the backlit LCD screen clearly displays the entered treatment time.
- 4 The durable slide-on comb attachment separates the hair for effective scalp penetration and maintains the proper distance from the skin to ensure consistent treatments.
- 5 The Phothera 100XL's treatment area is significantly larger than it's counterpart, the Phothera 100, making it an ideal solution for medium spot treatments.

Treatment Area	3.3" x 5.4"	8 cm x 14 cm/113 cm ²
Wand Size (Folded)	9" x 4.1" x 3"	23 cm x 10 cm x 8 cm
Wand Size (Open)	15" x 4.1" x 3"	38 cm x 10 cm x 8 cm
Wand Weight	1 lb	0.46 kg
Controller Size	7.25" x 2.75" x 4"	18.4 cm x 6.9 cm x 10 cm
Total Weight	3.25 lbs	1.47 kg

For further information, please contact Phothera:
216.831.0600
phothera.com | info@phothera.com