

# PhoThera<sup>PRO</sup> 4800 MAX



## PhoThera Pro 4800 MAX

Combining industry-leading technology with modern sophistication, the NeoLux has all the features a modern office is looking for. From multi-machine setup and customizable lamp configurations to advanced safety features, the PhoThera Pro 4800 MAX packs a powerful phototherapy punch in a small footprint.

PhoThera

# Pho<sup>o</sup>thera<sup>PRO</sup> 4800<sup>MAX</sup>



- 1 The open top and spacious interior helps to comfort claustrophobic patients and accommodate a wheel chair, while the stylish exterior will seamlessly blend into any office aesthetic.
- 2 The Pho<sup>o</sup>thera 4800 Pro Max offers ClearLink Controller capability via a 7" touch-screen controller, making treatment delivery quick, simple, safe, and effective.
- 3 48, high-intensity, six-foot lamps ensure treatment times are quick and effective.
  - Lamp Options Include: NBUVB, UVA, or a combination of the two
- 4 A removable platform elevates the patients to the desired height for effective full-body treatments.
- 5 UV-transmitting lamp shields keep dust and debris from accumulating and prevent patients from coming in contact with the lamps – and, they work in tandem with the lamp cooling fans to optimize the lamps operating temperature.
- 6 Patient-controlled fans offers a cooling breeze during treatments, ensuring patients stay comfortable and relaxed.
- 7 The Pho<sup>o</sup>thera Pro 4800 Max can be equipped with SmartTouch Technology and connectivity – offering remote connectivity, multi-machine setups, global and custom protocol creation, multi-language configuration, and more.

Height	84"	213 cm
Width (doors closed)	41"	104 cm
Width (doors open)	53"	135 cm
Depth (doors closed)	40"	102 cm
Depth (doors open)	50"	27 cm

For further information, please contact Pho<sup>o</sup>thera:  
216.831.0600  
pho<sup>o</sup>thera.com | clinical@pho<sup>o</sup>thera.com